


OPERATING INSTRUCTIONS


NOTE: Prior to Using your Hearing Aid, ensure that the unit is fully charged per the charging instructions on the reverse side.

Remove the Hearing Aid from the Case


1. Open the charging case.



2. Rotate the Hearing Aid vertically.




3. Remove the Hearing Aid from the case



4. Close the lid of the charging case.

Turning the Hearing Aid On and Off

• Use the Volume Adjustment Dial to turn the Hearing Aid On and Off.




1. To turn the power on, rotate the Volume Adjustment Dial clockwise to the desired volume level.

2. To turn the power off, rotate the Volume Adjustment Dial counter-clockwise until it stops.

Using the Hearing Aid

1. Turn the power on and ensure that the Ear Cap is securely fixed to the Hearing Aid.



2. While covering the microphone hole, place the Hearing Aid in your ear. **Not covering the microphone hole will result in a ringing sound.** If inserting the Hearing Aid into your right ear, the Volume Adjustment Dial should be up. If inserting the Hearing Aid into your left ear, the Volume Adjustment Dial should be down.

3. Use the Volume Adjustment Dial to adjust the volume to a comfortable level.

4. At the end of the intended use period or when a low battery chime warns you that the Hearing Aid needs to be recharged, refer to the charging instructions on the reverse side of this manual or the DVD for directions.

Troubleshooting

Headset will not charge:

- Are the two AAA batteries dead?
Try replacing the batteries or using the AC Adaptor to recharge the Hearing Aid.
- Is the Hearing Aid seated into the charging case correctly?
Ensure that the charging terminals are facing towards the LED light bar and that the Hearing Aid has been rotated down into the charging case correctly. View the enclosed DVD for a demonstration.
- Is the Hearing Aid on? (Note: Hearing Aid will not charge if in OFF Mode)
Rotate the Volume Adjustment Dial clockwise to ensure the hearing aid is on, then attempt to charge it again.

Can't hear anything:

- Is the power on?
Rotate the Volume Adjustment Dial clockwise to turn the power on.
- Is the Ear Cap blocked with wax or debris?
Clean the Ear Cap using the provided cleaning brush. Do not wash the ear cap in water or use any industrial cleaners or solvents on the ear cap.
- Is the rechargeable battery dead?
Place the Hearing Aid in the charging case for 8-12 hours and ensure it is charging properly.

Warnings and Notices

Warning to Hearing Aid Dispensers

A hearing aid dispenser should advise a prospective hearing aid user to consult promptly with a licensed physician (preferably an ear specialist) before dispensing a hearing aid if the hearing aid dispenser determines through inquiry, actual observation, or review of any other available information concerning the prospective user, that the prospective user has any of the following conditions:

- Visible congenital or traumatic deformity of the ear.
- History of active drainage from the ear within the previous 90 days.
- History of sudden or rapidly progressive hearing loss within the previous 90 days.
- Acute or chronic dizziness.
- Unilateral hearing loss of sudden or recent onset within the previous 90 days.
- Audiometric air-bone gap equal to or greater than 15 decibels at 500 hertz (Hz), 1,000 Hz, and 2,000 Hz.
- Visible evidence of significant cerumen accumulation or a foreign body in the ear canal.
- Pain or discomfort in the ear

Special Care should be exercised in selecting and fitting a hearing aid whose maximum sound pressure level exceeds 132 decibels because there may be risk of impairing the remaining hearing of the hearing aid user. (This provision is required only for those hearing aids with a maximum sound pressure capability greater than 132 decibels (dB).) This Hearing Aid does not exceed 125 decibels.

Important Notice for Prospective Hearing Aid Users

Good health practice requires that a person with a hearing loss have a medical evaluation by a licensed physician (preferably a physician who specializes in diseases of the ear) before purchasing a hearing aid. Licensed physicians who specialize in diseases of the ear are often referred to as otolaryngologists, otologists or otorhinolaryngologists. The purpose of medical evaluation is to assure that all medically treatable conditions that may affect hearing are identified and treated before the hearing aid is purchased.

Following the medical evaluation, the physician will give you a written statement that states that your hearing loss has been medically evaluated and that you may be considered a candidate for a hearing aid. The physician will refer you to an audiologist or a hearing aid dispenser, as appropriate, for a hearing aid evaluation.

The audiologist or hearing aid dispenser will conduct a hearing aid evaluation to assess your ability to hear with and without a hearing aid. The hearing aid evaluation will enable the audiologist or dispenser to select and fit a hearing aid to your individual needs.

If you have reservations about your ability to adapt to amplification, you should inquire about the availability of a trial-rental or purchase-option program. Many hearing aid dispensers now offer programs that permit you to wear a hearing aid for a period of time for a nominal fee after which you may decide if you want to purchase the hearing aid.

Federal law restricts the sale of hearing aids to those individuals who have obtained a medical evaluation from a licensed physician. Federal law permits a fully informed adult to sign a waiver statement declining the medical evaluation for religious or personal beliefs that preclude consultation with a physician. The exercise of such a waiver is not in your best health interest and its use is strongly discouraged.

Children With Hearing Loss

In addition to seeing a physician for a medical evaluation, a child with a hearing loss should be directed to an audiologist for evaluation and rehabilitation since hearing loss may cause problems in language development and the educational and social growth of a child.

NOTICE

- * This product is not intended for use by anyone under the age of 18.
- * For mild to moderate hearing loss.
- * Hearing Aids are not designed to restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions.
- * In most cases, infrequent use of a hearing aid does not permit a user to attain full benefit from it.
- * The use of a hearing aid is only part of hearing habilitation and may need to be supplemented by auditory training and instruction in lip reading.
- * Do not attempt to use the hearing aid in a method not instructed within this manual.
- * Keep this product away from infants, children and pets.
- * Remove your hearing aid prior to sleeping.
- * Do not use this product while operating a mobile phone if it is uncomfortable for you to do so. Radio wave disturbances may occur.
- * Remove your hearing aid prior to vigorous exercise to avoid injury or damage to yourself or to the hearing aid.
- * Remove your hearing aid when exposure to moisture is possible (showering, washing your face, using hair spray, etc.)
- * Please consult local authorities and regulations regarding the disposal of your hearing aid and the internal rechargeable battery.

SPECIFICATIONS

General Product Description:	Rechargeable Hearing Aid
Title of Commodity / Item Number:	Hearing Aid / UP-6S86
Standard Frequency:	1600Hz
Frequency Range:	between 200 and 6000 Hz
Max.output sound level gain:	29 ± 5dB 0~29dB
Highest sound pressure level for input SSPL 90dB:	118 ± 5dB below 125dB
Equivalent input noise level:	below 29dB
Reference Test Gain:	above 30dB
Average Saturation Output (HF-Average SSPL 90):	113 ± 5db
Saturation Output Curve (SSPL 90 Curve):	118 ± 5db
Frequency Response Curve:	between 200 and 6000 Hz
Total Harmonic distortion value :	1600Hz: below 3%;800Hz:below 4%;500Hz:below 4%
Battery Current Drain :	below 0.8mA / 2.5V
Battery Type / life expectancy:	button size rechargeable battery. continuous use: 6~10 hours per charge. 500+ recharges
Operating temperature / humidity :	+10°C ~ +40°C (50F-104F),30% ~ 85% humidity
Storage temperature / humidity :	-10°C ~ +60°C (50F-140F),30% ~ 85% RH
Size (L/W/H):	<Hearing Aid> : 10.5 x 14.4 x 24.0 mm / 0.41 x 0.57 x 0.94 in <Charging Case>: 96.0 x 73.0 x 30.0 mm / 3.78 x 2.87 x 1.18 in
Weight :	<Hearing Aid>2.4g <Charging Case>67g
Material :	<Hearing Aid> ABS <Charging Case>ABS <Ear Cap> Silicone <AC Adapter>ABS